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PEG, percutaneous endoscopic gastrostomy

References:

- Elisa T, Adrian T, Yu-Ling H. Gastrostomy tube feeding and management in the community. Singapore Fam Physician 2020; 46(6): 27-33.
- Cardiovascular and Interventional Radiological Society of Europe. Gastrostomy [Internet]. [cited 2021 Jul 27]. Available from: <https://www.cirse.org/patients/ir-procedures/gastrostomy/>.
- Warriner L, Spruce P. Managing overgranulation tissue around gastrostomy sites. Br J Nurs. 2012; 21(Sup5):S14-24.
- Best C. The correct positioning and role of an external fixation device on a PEG. Nursing times. 2004; 100(18):50-1.
- A Clinician's Guide to Caring for people with gastrostomy tubes and devices. AGENCY FOR CLINICAL INNOVATION (ACI) guidelines. Available from: https://www.aci.health.nsw.gov.au/_data/assets/pdf_file/0017/251063/gastrostomy_guide-web.pdf.
- Malhi H, Thompson R. PEG tubes: dealing with complications. Nurs Times. 2014 Nov 1;110(45):18-21.
- Rahnamaiazar AA, Naghshizadian R, Kurtz A, Farkas DT. Percutaneous endoscopic gastrostomy: indications, complications and management. World J Gastroenterol. 2014;20(24):7739-7751.
- Blumenstein I, Shastri YM, Stein J. Gastroenteric tube feeding: techniques, problems and solutions. World J Gastroenterol. 2014;20(26):8505-8524.
- Avanos. MIC-KEY* Feeding tube kits [Internet]. [cited 2021 Jul 27]. Available from: <https://avanosmedicaldevices.com/digestive-health/enteral-feeding/mic-key-feeding-tube-kits/>.
- Product data sheet. HC205-01_DH Legacy_MIC-KEY_threefold_UK_LR.
- Product data sheet. HC113-04-UK_MIC&MIC-KEY_EternalFeeding_ProductCatalogue_2020.pdf.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2018 AVNOS. All rights reserved. COPY-05620