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References:

1. Shah J, Shahidullah A, Richards S. Reducing the unintended dislodgement of gastrostomy tubes in a long-term acute care hospital: A QA/QI pilot study. Gastroenterol Res. 2018;11(5):369-373.
2. Sealock RJ, Munot K. Common gastrostomy feeding tube complications and troubleshooting. Clinical Gastroenterology and Hepatology. 2018;16(12):1864-9. 3. Feil M. Dislodged gastrostomy tubes: preventing a potentially fatal complication. Reviews Analyses. 2017;14(1):9-16. 4. Tips To Avoid G-Tube Dislodgement. https://tubefed.com/newsletter/tips-to-avoid-g-tube-dislodgement/. Accessed on 27.08.2021. 5. Rowat A. Enteral tube feeding for dysphagic stroke patients. British Journal of Nursing. 2015 Feb 12;24(3):138-45. 6. Schrag SP, Sharma R, Jaik NP, Seamon MJ, Lukaszczuk JJ, Martin ND, Hoey BA, Stawicki SP. Complications related to percutaneous endoscopic gastrostomy (PEG) tubes. A comprehensive clinical review. Journal of Gastrointestinal and Liver Diseases. 2007 Dec 1;16(4):407. 7. Rahnama-Azar AA, Rahnamaiazar AA, Naghshizadian R, Kurtz A, Farkas DT. Percutaneous endoscopic gastrostomy: indications, technique, complications and management. World J Gastroenterol. 2014;20(24):7739-7751. 8. Boullata JI, Carrera AL, Harvey L, Escuro AA, Hudson L, Mays A, McGinnis C, Wessel JJ, Bajpai S, Beebe ML, Kinn TJ. ASPEN safe practices for enteral nutrition therapy. Journal of Parenteral and Enteral Nutrition. 2017 Jan;41(1):15-103. 9. MIC* PEG sell sheet copy-04844. 10. AVANOS catalogue 2020 ANZ-DH-New. 11. MIC-KEY* G Feeding Tubes.Product data sheet. HC205-01_DH Legacy_MIC-KEY_threefold_UK_LR. 12. Product data sheet, MIC* and MIC-KEY* enteral feeding product catalogue 2020.